

MTB/NTM Screening & Drug-Resistant TB Detection

using Real-time PCR

CE-IVD Marked



- Ⓐ Anyplex™ MTB/NTM Real-time Detection**
 - Simultaneous detection of MTB and NTM
- Ⓐ Anyplex™ II MTB/MDR/XDR Detection**
 - Simultaneous screening for MTB infection, 25 MDR mutations and 13 XDR mutations
- Ⓐ Anyplex™ II MTB/MDR Detection**
 - Simultaneous detection of MTB and 25 mutations associated with MDR-TB
 - Rifampicin-resistance (18 mutations)
 - Isoniazid-resistance (7 mutations)
- Ⓐ Anyplex™ II MTB/XDR Detection**
 - Simultaneous detection of MTB and 13 mutations associated with XDR-TB
 - Fluoroquinolone-resistance (7 mutations)
 - Injectable drug-resistance (6 mutations)



HIGH SENSITIVITY & SPECIFICITY

Multiplex real-time PCR with high sensitivity and specificity by utilization of DPO™ and TOCE™ technologies



Cost-effective assays for MTB, NTM and drug-resistant TB

Anyplex™ MTB/NTM Real-time Detection is an efficient method to screen TB quickly in the first instance.

Anyplex™ II MTB/MDR, Anyplex™ II MTB/XDR and Anyplex™ II MTB/MDR/XDR Detection are able to identify TB and drug-resistant TB simultaneously and rapidly, which allows cost-saving test. An ideal testing algorithm is to screen patient samples with Anyplex™ MTB/NTM Real-time Detection followed by three assays. It allows reductions in cost, time, and labor for additional sample collection and nucleic acid extraction.

Complete Solution for TB Screening and Drug-Resistant TB Detection


- Enable to distinguish or identify TB and drug-resistant TB within clinically meaningful timeframe
- Eliminate need for additional sample collection and nucleic acid extraction to test drug-resistant TB when coupled with MTB/NTM screening
- Provide prompt diagnosis and appropriate treatment guideline for TB control

Features

1. Available for various specimen types

- Sputum
- Cultured cell
- Bronchial washing
- Fresh tissue

2. Reportable within one working day of sample receipt



* Turnaround time : nucleic acid extraction to result

3. Convenient result analysis and interpretation through Seegene Viewer

- Interface specialized for multiplex testing
- Interlocked with LIS
- Easy readout for comprehensive information

Workflow

- Collection of sample from suspected TB patient



Nucleic acid extraction & PCR Setup

CFX96™ Dx



Real-time PCR

Seegene Viewer



Automated data analysis

MTB / MDR Detection

Simultaneous detection of *Mycobacterium tuberculosis* (MTB) and 25 mutations associated with MDR-TB



Analytes

- **MTB**
- **Internal Control (IC)**
- **Multi-Drug Resistance (MDR)**
 - Isoniazid-resistance (7 mutations)
 - Rifampicin-resistance (18 mutations)


1. Broad coverage of MDR-TB point mutations

		Target mutations of MDR-TB detection				
Drug resistance	Related gene	Target mutations				
RIF-R* (18 mutations)	<i>rpoB</i>	L511P (CTG ⇒ CCG)	3 a.a. deletion in 513~516	S522Q (TCG ⇒ CAG)	H526N (CAC ⇒ AAC)	S531W (TCG ⇒ TGG) L533P (CTG ⇒ CCG)
		Q513K (CAA ⇒ AAA)	D516V (GAC ⇒ GTC)	H526C (CAC ⇒ TGC)	H526R (CAC ⇒ CGC)	
		Q513L (CAA ⇒ CTA)	D516Y (GAC ⇒ TAC)	H526D (CAC ⇒ GAC)	H526Y (CAC ⇒ TAC)	
		Q513P (CAA ⇒ CCA)	S522L (TCG ⇒ TTG)	H526L (CAC ⇒ CTC)	S531L (TCG ⇒ TTG)	
INH-R (7 mutations)	<i>katG</i>	S315I (AGC ⇒ ATC)	S315N (AGC ⇒ AAC)	S315T (AGC ⇒ ACC)	S315T (AGC ⇒ ACA)	
	<i>inhA</i> promoter	-15 (C ⇒ T)	-8 (T ⇒ A)	-8 (T ⇒ C)		

* It is possible to detect 9 additional RIF resistance mutations, which have the same codon site.

2. Anyplex™ II MTB/MDR Detection provides more information for appropriate treatment¹⁻³⁾

- Schematic diagram of test result



- Clinical implication according to melting temperature (INH-R)

Interpretation	Melting temp. (Tm)	Corresponding gene / mutations	Clinical implication (general aspect)
INH-R	1 Low Tm	4 mutations in <i>katG</i>	High-level INH resistance
	2 High Tm	3 mutations in <i>inhA</i> promoter	Low-level INH resistance


Result / Seegene Viewer

References

1. Guo H. *et al.*, J. Med. Microbiology (2006) 55:1527-31
2. Johnson R. *et al.*, Curr Issues Mol Biol. (2006) 8:97-111
3. Ando H. *et al.*, Antimicrob Agents Chemother. (2010) 54:1793-9

MTB / XDR Detection

Simultaneous detection of *Mycobacterium tuberculosis* (MTB) and 13 mutations associated with XDR-TB



Analyses


- MTB
- Internal Control (IC)
- Extensively Drug Resistance (XDR)
 - Fluoroquinolone-resistance (7 mutations)
 - Injectable drug-resistance (6 mutations)

1. Broad coverage of XDR-TB point mutations

Target mutations of XDR-TB detection						
Drug resistance	Related gene	Target mutations				
FQ-R (7 mutations)	<i>gyrA</i>	A90V (GCG ⇒ GTG) S91P (TCG ⇒ CCG)	D94A (GAC ⇒ GCC) D94G (GAC ⇒ GGC)	D94H (GAC ⇒ CAC)	D94N (GAC ⇒ AAC)	D94Y (GAC ⇒ TAC)
Injectable drug-R (6 mutations)	<i>rrs</i>	1401 (A ⇒ G)	1402 (C ⇒ T)	1484 (G ⇒ T)		
	<i>eis</i> promoter	-37 (G ⇒ T)	-14 (C ⇒ T)	-10 (G ⇒ A)		

2. Anyplex™ II MTB/XDR Detection provides more information for appropriate treatment⁴⁻⁶⁾

- Schematic diagram of test result




- Clinical implication according to melting temperature (Injectable drug-R)

Interpretation	Melting temp. (Tm)	Corresponding gene / mutations	Clinical implication (general aspect)
Injectable drug-R	1 Low Tm	3 mutations in <i>eis</i> promoter	Low-level KAN resistance
	2 Middle Tm	2 mutations in <i>rrs</i> (1401G/1484T)	High-level KAN/AMI/CAP resistance
	3 High Tm	1 mutation in <i>rrs</i> (1402T)	Low-level KAN resistance High-level CAP resistance

KAN : Kanamycin, AMI : Amikacin, CAP : Capreomycin

Result / Seegene Viewer



References

- Johnson R. et al, Curr Issues Mol Biol. (2006) 8:97-111
- Zaubrecher MA. et al, Proc Natl Acad Sci USA (2009) 106:20004-9
- Gikalo MB et al, J Antimicrob Chemother. (2012) 67:2107-9



REAL TIME CE-IVD Marked



MTB / NTM Real-time Detection

Simultaneous detection of *Mycobacterium tuberculosis* (MTB) and non-tuberculosis mycobacteria (NTM)

Analyses

- MTB
- Mycobacteria
- Internal Control (IC)

Result / Seegene Viewer



REAL TIME CE-IVD Marked


MTB / MDR / XDR Detection

Simultaneous detection of MTB, 25 mutations associated with multi-drug-resistant tuberculosis (MDR-TB) and 13 mutations associated with extensively drug-resistant tuberculosis (XDR-TB)

Analyses

- MTB
- Internal Control (IC)
- Multi-Drug Resistance
 - Isoniazid-resistance (7 mutations)
 - Rifampicin-resistance (18 mutations)
- Extensively Drug Resistance
 - Fluoroquinolone-resistance (7 mutations)
 - Injectable drug-resistance (6 mutations)

Result / Seegene Viewer



Benefits of Seegene's TB Solution



For Patients

- Be treated promptly
- Get personalized treatment



For Physician Clinics

- Give faster answers to patients
- Aid in selecting appropriate therapy (first-line or second-line anti-TB drugs)



For Technicians / Labs

- Reduce hands-on time with health-threatening specimen
- Provide comprehensive tests for MTB/NTM and MTB/MDR/XDR on one platform



For TB Controller

- Get more results (MTB, NTM, MDR and XDR) to deliver timely and targeted therapy
- Control secondary transmission of MTB and drug-TB resistance in the community



• Ordering information

Not Available for Sale in the United States

Product	Package Volume	Cat. No.
Anyplex™ MTB/NTMe Real-time Detection	50 rxns	TB7202Y
	100 rxns*	TB7202X
Anyplex™ MTB/NTM Real-time Detection (V2.0)	100 rxns	TB7200X
Anyplex™ II MTB/MDR/XDR Detection	50 rxns	TB7500Y
Anyplex™ II MTB/MDR Detection	50 rxns	TB7301Y
Anyplex™ II MTB/XDR Detection	50 rxns	TB7302Y
Instrument	Type	Cat. No.
CFX96™ Dx	Real-time PCR _ Optical Reaction Module	1845097-IVD
	Real-time PCR _ Thermal Cycler	1841000-IVD

*For use with NIMBUS IVD & STARlet IVD only



Taewon Bldg. 91 Ogeum-ro, Songpa-gu, Seoul 05548, Republic of Korea / Tel : +82-2-2240-4000 / Fax : +82-2-2240-4040 / E-mail : info@seegene.com

www.seegene.com

BRAZIL

Belo Horizonte, Brazil
Tel : +55-31-25153003
E-mail : contato@seegenebrazil.com.br

CANADA

Toronto, Canada
Tel : +1-800-964-5680
E-mail : canada@seegene.com

GERMANY

Düsseldorf, Germany
Tel : +49-211-9943-4260
E-mail : sgg@seegene.com

MEXICO

México city, México
Tel : + 52 (55)-8848-9646
E-mail : mexico@seegene.com

MIDDLE EAST

Dubai, UAE
Tel : +971-4-558-7110
E-mail : sgme@seegene.com

USA

California, USA
Tel : +1-925-448-8172
E-mail : usa@seegene.com