



ALLIANCE
BIO EXPERTISE
Connected Microbiology



AIRWEL[®]

Bio-collector





AIRWEL®


*Increase your
productivity
quality
traceability*

**ADAPTABLE TO
YOUR NEEDS**

PETRI DISH:
90 mm
55 mm

AIRFLOW:
100L / min
200L / min

FULL TRACEABILITY
GMP / GLP – 21 CFR PART 11



EXCLUSIVE OPTION

IN SITU AIRFLOW CHECK

Check your airflow
whenever you like !

Don't wait the 1 year
calibration services to check

Avoid wrong collection due
to wrong airflow !

SAFETY MANAGEMENT

Cascade password

3 level:

Operator – Manager – Service



FULLY HERMETIC
ip65 dust and water protection

WIRELESS BATTERY CHARGER

ANTIBACTERIAL PAINT

SERVICES

MANUFACTURER or DEALER **CALIBRATION**
QUALIFICATION: **IQ/OQ** protocol

CONNECTED
MICROBIOLOGY





SPECIFICATIONS

AIRWEL®



COMPLIANCE	AIRWEL	AIRWEL+
IP65		✓
ISO14698 / ISO17025	✓	✓
USP797 / USP116	✓	✓
EU GMP / GLP	✓	✓
21 CFR PART 11	✓	✓
TRACEABILITY		
Bluetooth module	✓	✓
PC software for data management	✓	✓
Android App (PC / Smartphone / tablet)	✓	✓
CALIBRATION		
IN SITU AIRFLOW CHECK (Option)		✓
A.B.E or dealer calibration	✓	✓
BENEFITS		
Cascade password		✓
Wireless Induction charger		✓
Battery autonomy	More than 10 hours	
Type of dishes	90 mm or 55 mm	
Type of heads	Stainless steel or Disposable polyester	
Airflow	100L or 200L / min	
Transport case	Light	Rigid
Multi-mode aspirating cycles	Uninterrupted	Fragmented
Antigerm paint	✓	✓
DIMENSIONS		
W x H x D	130 x 330 x 150 mm	
Weight	1.6 kg	
Power supply / Electric power	120-230 V – 50-60HZ / 120W	
Made in EU	✓	✓



ALLIANCE BIO EXPERTISE

Empowered from our people skills and experiences, we, ALLIANCE BIO EXPERTISE, together with our customers and partners decided to manufacture Microbiology automated equipment around essential values :
Customer Focus, Service, Loyalty, Innovations

**Ask your local partner
for a personal demonstration!**

ZI de courbouton , Le Tremplin – 35480 GUIPRY - FRANCE
E-mail : contact@abioexpertise.com
Phone : +33(0)2.40.51.79.53
alliance-bio-expertise.com

Follow us on 

