




⊕ THE FIRST COMPACT & BENCH TOP FULLY AUTOMATED INSTRUMENT FOR  
**IMMUNOFIXATION** ON AGAROSE GEL FROM PRIMARY TEST TUBE




Discover more  
on the EASYFIX

⊕ The new Easyfix Interlab G26 performs all electrophoresis  
assays from **primary tubes** in full automation from  
the conventional serum proteins to immunofixation

*EASYFIX Interlab G26 one step ahead!!!*



Automated sample  
dilution



Automated sample  
application



Automated antisera  
application



Automated  
blotting



Automated  
Gel reader



## ALL IMMUNOFIXATION PHASES IN **TOTAL AUTOMATION** WITHOUT ANY OPERATOR INTERVENTION

- BARCODE READER: positive Sample ID from Primary test tube
- GEL ID: UNIQUE IDENTIFIER for sample tracking
- ASSAY ID: Assay specific identifier communicating procedure directly to the instrument
- CONTINUOUS SAMPLE PROCESSING: Primary tube sample rack allowing continuous specimen addition for non-stop processing.
- FLEXIBILITY: Different or same assays can be processed simultaneously (SPE, IFE)
- FAST: Multiple gel configurations to handle all workloads. Up to 12 IFE and 117 SPE per hour
- EXTENSIVE MENU: 15 available procedures, from SPE to IEF of the CSF.
- COMPACT: Just 85 x 80 x 53 cm

**Interlab Srl**

Via Rina Monti 26, 00155 Rome, Italy | Tel: +39-06-22754350 | Fax: +39-06-22754534 | e-mail: [info@interlab-srl.com](mailto:info@interlab-srl.com)

<http://www.interlab-srl.com>