

minicap sebia
flex piercing

The capillary
electrophoresis
based on a solid
experience



minicap sebia

flex piercing

More flexibility for your laboratory

Multi-parameter instrument

Serum, whole blood capped tubes

Large test menu

Proteins, Immunotyping, CDT, Hemoglobin and HbA1c

Fully automated

Cap piercing: maximizes productivity and reduces biohazard handlings

Positive ID: full traceability from the primary tube to final result

Analytical superiority

Assay parameters optimized for each technique

High resolution

Clear-cut and precise separation

Throughput / 2 capillaries

20 tests/h (Protein)

8 tests/h (HbA1c and Hb)

10 tests/h (CDT)

Possible adaptation of the MINICAP into MINICAP Flex Piercing

Adaptation kit available on request (Cat.Nb 1240)



sebia

Parc Technologique Léonard de Vinci
CP 8010 Lisses - 91008 EVRY Cedex
France

sebia Benelux SCS / Comm.V.

Jan Olieslagerslaan, 41

1800 Vilvoorde

Belgique / België

Tel.: +32 (0)2 702 64 64

Fax: +32 (0)2 702 64 60

e-mail: sebia.benelux@sebia.be

sebia Brasil Ltda

Edifício Baker Office Tower,

Rua Barão do Triunfo, 73, Conjunto 51

Bairro Brooklin Paulista, CEP 04602-000

Brasil

Tel.: +55 11 3849 01 48

Fax: +55 11 3841 98 16

e-mail: sebia@sebia.com.br

sebia GmbH

Münsterfeldallee, 6

36041 Fulda

Deutschland

Tel.: +49 (0)661 3 30 81

Fax: +49 (0)661 3 18 81

e-mail: sebia@sebia.de

sebia Hispania S.A.

C / Sicilia, n° 394

08025 Barcelona

España

Tel.: +34 93 208 15 52

Fax: +34 93 458 55 86

e-mail: sebia@sebia.es

sebia Inc.

400-1705 Corporate Drive

Norcross, GA 30093

U.S.A.

Tel.: +1 770 446 - 3707

Fax: +1 770 446 - 8511

e-mail: info@sebia-usa.com

sebia Italia S.r.l.

Via Antonio Meucci, 15/a

50012 Bagno a Ripoli (FI)

Italia

Tel.: +39 055 24851

Fax: +39 055 2485400

e-mail: info@sebia.it

sebia UK Ltd

River Court, Meadows Business Park

Station Approach, Blackwater

Camberley, Surrey, GU17 9AB

United Kingdom

Tel.: +44 (0)1276 600636

Fax: +44 (0)1276 38827

e-mail: sales@sebia.co.uk

sebia Shanghai Rep. office

Cross Tower, Room 2306-07,

318 Fuzhou Road

Shanghai 200001

China

Tel.: +86(21) 6350 1655

Fax: +86(21) 6361 2011

e-mail: sebia@sebia.cn

