

PETRI DISH WRAPPER

Pack your Petri dishes

MINIMA®

Order it with your Distri wel®



Adequate film
optimal solution to avoid condensation



..... Packs your petri dishes in packs of 10
(one column of 10 or two columns of 5)

Programs
10 predefined



Wrapping machine Minima - Packs your petri dishes



Wrapping film (delivred with)	Polyolefin rolls delivered by 6 rolls of 1250m – 15µm – Ø 200mm
Power supply	200/230V - 1.2 KW
Dimensions with open lid	930 x 550 x 580 mm
Dimensions with closed lid	930 x 550 x 370 mm
Weight (net/gross)	36/48 kg

