

**RIEGER**

Industrievertretungen GmbH

www.rieger-iv.at

Wir sind für Sie da:

E-Mail: office@rieger-iv.at

Telefon: +43(1) 728 00 52

Fax: +43(1) 728 69 16

Rustenschacher Allee 10

1020 Wien, Austria

MIXWEL®

400 mL lab blenders





MIXWEL®

Increase your
productivity
quality
traceability

PERFORMANCES



Silent <48db



28 to 70 kg
pressure
per paddle

CONNECTED
MICROBIOLOGY



PATENTED pendular blending
Delivers up to
70 kg pressure /paddle.
Optimum homogenization !

Extreme Compactness

The most compact blender
on the market (30 dm3)

Space optimization of your bench

Easily transportable

Comfort of cleanability

Quick removal of the door

Removable paddles

Waste tank integrated
into the door



Stainless
steel chamber (AISI 304)



**Lightest of
the market !**
14 kg

VIBRATION FREE
Extreme stability
while blending



SPECIFICATIONS

MIXWEL®



SPECIFICATIONS

Pendular blending power	Up to 70 kg /paddle
Capacity	50 mL to 400 mL
Noise level	< 48 dbA
Speed	380 rpm

ADVANTAGES

Vibration free	✓
Removable paddles	✓
Stainless steel paddles	✓
Adjustable paddles to the matrix	3 positions
Waste tank integrated	✓

DIMENSIONS & ELECTRICAL SUPPLY

Compacity W x D x H	220 x 410 x 340 mm
Weight	14 kg
Electrical supply	230 V – 50 Hz
Power	90 W

THEY HAVE ADOPTED IT

"I didn't know we could blend much harder
in a smooth way" ... " We blend young elastic
cheese ... and the result is perfect, similar
to milk aspect with no solids at all."

Lab user, INRA STLO RENNES

"As the microorganisms' extraction is much
better with the MIXWEL, we have been able to
reduce the blending time" ... "never drilled a bag" ...
"a quiet blender, small, with no vibrations"

Mrs Ali Mandjee, RAEMA / ASA PARIS Region



Empowered from our people's skills and experiences, we, ALLIANCE BIO EXPERTISE, together with our customers and partners decided to manufacture Microbiology automated equipment around essential values:
Customer Focus, Service, Loyalty, Innovations

**Ask your local partner
for a personal demonstration!**

ZI de courbouton , Le Tremplin – 35480 GUIPRY - FRANCE
E-mail : contact@abioexpertise.com
Phone : +33(0)2.40.51.79.53
alliance-bio-expertise.com