

POLYWEL 3D[®]

AUTOMATIC TUBES AND BOTTLES FILLER

Versatility, accuracy and speed




AUTOMATED TUBES AND BOTTLES FILLER

Versatility, accuracy and speed

POLYWEL 3D[®]



1000 tubes per hour / per run*



8 racks of 120 tubes or 4 racks of 250 tubes*

3D driven

Flexible and adaptable to your different containers

Distribution nozzles :

Ø 1mm - Ø 3mm - Ø 6mm

ACCURACY → SPEED

Autoclavable racks

To be loaded directly from its original package

End of the cycle light indicator

Up to 20 programs, quick to start

- Accuracy from 5ml : ± 0.1 ml
- Calibration at the beginning of the run
- Calibration during the run

Easy cleaning

stainless steel surface






* Ø 18mm tubes (9 ml)



Delivered with the POLYWEL 3D : User Manual, power chord, Polypump 3D (POLW2000), Set of tubing 3.2 mm with dispensing nozzle (POLW2008), Dispensing nozzle holder (POLW2013), rack 250.

TECHNICAL DATA SHEET

Distribution volume	From 1 to 9999.9ml
Tubes per hour (Ø 18mm)	1000
Tubes per rack (Ø 18mm)	Large rack : 250 / Small rack : 120
Programs	20
Tubing diameters	1.6 mm / 3.2 mm / 6.4 mm
Adjustable speed	✓
Calibration	Start and during the cycle
Maximum size of containers	Height : between 115 and 250mm Width : up to 190mm Depth : up to 390mm
Dimensions	1130 x 600 x 610 mm
Power supply	120 – 240V (±10%) 50 – 60 Hz / 160W
Warranty	1 year
Country of origin	European Union - France



POLYWEL 3D®

Distribution of media, diluent or other liquids in tubes, bottles, Petri dishes, square dishes ...



ACCESSORIES



Delivred with :

Polywel 3D (POLW1001) - Polypump 3D (POLW2000) - Set of tubing 3.2mm with distribution nozzle (POLW2008) - Rack 250 (POLW2002) - Connection cable - USB key with user manual and quality report

Accessories :

- POLW2000 : Polypump 3D
- POLW2001 : Rack 120
- POLW2002 : Rack 250
- POLW2003 : Angled loading holder rack holder
- POLW2004 : Set of tubing 1.6 mm without distribution nozzle
- POLW2005 : Set of tubing 3.2 mm without distribution nozzle
- POLW2006 : Set of tubing 6.4 mm without distribution nozzle
- POLW2007 : Set of tubing 1.6 mm with distribution nozzle
- POLW2008 : Set of tubing 3.2 mm with distribution nozzle
- POLW2009 : Set of tubing 6.4 mm with distribution nozzle
- POLW2010 : Dispensing nozzle 1.6 mm
- POLW2011 : Dispensing nozzle 3.2 mm
- POLW2012 : Dispensing nozzle 6.4 mm
- POLW2013 : Dispensing nozzle holder
- POLW2014 : Accuracy pump head



**20 YEARS EXPERIENCE
IN MICROBIOLOGY EQUIPMENT DESIGN**

**A FULLY COMMITTED TEAM
OF SERVICE ENGINEERS AND MICROBIOLOGISTS**



**MORE THAN 70 INTERNATIONAL DISTRIBUTORS
TRAINED TO THE LATEST A.B.E. TECHNOLOGIES**

**A CONNECTED EQUIPMENT GENERATION
FOR TRACEABILITY AND PRODUCTIVITY**



FOOD SAMPLE PREPARATION


Blenders



Diluters




Tubes Fillers Pumps




AUTOMATED MEDIA PREPARATION

Media Preparators



Petri Fillers



SAMPLE COLLECTION

Air Samplers



Don Whitley Scientific Limited

Victoria Works, Victoria Street, Bingley, BD16 2NH, UK
t: +44 (0)1274 595728 e: sales@dwscientific.co.uk



**don whitley
scientific**
www.dwscientific.com