


WELD MC GRIFOLS

Desk-top multiple tube sealer with 4 automatic and programmable sealers.

- Maximum efficiency.
- Fast: sealing time 7 seconds.
- High safety, it incorporates anti-spark and overheat protection systems.



WELD T GRIFOLS

Automatic desk-top tube sealer for PVC medical grade tubing.

- Designed to withstand intensive use.
- Automatic sealer head.
- Automatically regulates the sealing time and pressure.



WELD B GRIFOLS

Portable tube sealer with battery for PVC medical grade tubing.

- Light-weight and easy to use.
- Ideal for locations that require high mobility.
- Up to 500 seals.



TECHNICAL CHARACTERISTICS

	Weld MC GRIFOLS	Weld T GRIFOLS	Weld B GRIFOLS
Power source	115/230VAC 50/60HZ	115/230VAC 50/60HZ	Internal battery 110/240VAC 47/63HZ
Consumption	550VA	420VA	700mA
Maximum RF power	80W	80W	80W
Dimensions (cm)	30Wx12.7Hx20.8D	20Wx12Hx22D	10.7Wx4.2Hx16.1D
Net weight	6KG	4KG	1KG